


Validation of Suprachoroidal Injection Training Program with a Synthetic Eye Model

Cherry Wan¹, Nathan Fisher¹, Rafael Andino¹
¹Clearside Biomedical, Alpharetta, GA, USA

PURPOSE

To demonstrate the effectiveness of a proprietary training program for suprachoroidal injection with the SCS Microinjector®, the device constituent of an FDA approved combination product, XIPERE®

Figure 1: SCS Microinjector Syringe with 900 µm Microneedle



BACKGROUND

- Over 15,000 suprachoroidal injections have been performed with the SCS Microinjector, in over 15 clinical trials (Table 1) and commercial use. It is the only drug delivery device constituent approved for suprachoroidal injection in a combination product by the FDA.
- As a new route of administration, a robust and standardized training procedure was developed to ensure consistency across multiple clinical trials, indications, products, and countries.
- A two-tiered proprietary training program was developed: Training the Trainers (TtT) and physician training (PT).
- To date, there are 62 global certified trainers and thousands of trained retinal specialists.
- Training material, including a custom-developed model eye and training content (Figure 2), has been developed to provide an effective training experience.

A new Category 1 CPT Code has contributed to expanded commercial use and adoption of suprachoroidal injections with the SCS Microinjector® with XIPERE® (Bausch + Lomb).

The suprachoroidal injection training program with the SCS Microinjector® is **THE ONLY CLINICALLY VALIDATED** training program and over 15,000 injections have been completed to date.

METHODS

- Three (3) surveys have been completed to evaluate user understanding of the procedure after training.
- Early Adopter Feedback: twelve (12) retinal/uveitis specialists with 10+ suprachoroidal injections of XIPERE completed a 37-question survey.¹
- A total of 33 people were trained and surveyed to evaluate the adequacy of the training program.²
- User survey for injecting physicians and observers from one of the Clearside Phase 3 studies for uveitis.³
- Critical tasks related to the procedure were determined with a User-Related Risk Analysis (URRA).⁴
- Effectiveness of the training was confirmed with device and procedural complaint data during multiple clinical trials with a total of 444 injections.
- Training model effectiveness impact was ranked from level 1 to 5, from trainee satisfaction to impact of procedure at a systemic health care level (Table 2).^{5,6}

Table 2: Translational outcome of simulation-based learning

Levels	Parameter	Definition	Example
Level 1	Internal acceptability	The trainee's satisfaction with using the simulator	Favorable responses from feedback forms or post-training survey questionnaires
Level 2	Contained effects	Changes in performance in the simulation context	Changes in performance in the simulation context
Level 3	Downstream effects	Behavioral change in the clinical context	Transfer of knowledge/skills to clinical practice
Level 4	Target effects	Direct changes to patient outcomes	Reduce rates of complications
Level 5	Collateral effects	Changes on a wider, systemic levels	Cost saving; skill retention


Disclosure: All authors (CW, NF and RA) are employees of Clearside Biomedical, Inc.



Presentation
5819-B0415



Figure 2: Proprietary training program: a custom training eye model with training material content (Included in figure is an excerpt of the critical tasks associated with the injection procedure)



RESULTS

- Early adopter survey showed **training improved physician comfort in injection procedure.**¹
- Over 64% of trainees believe **no mandated retraining was necessary**, with the rest recommending annual re-training or refresher.²
- Clinical trial questionnaire showed **84% of physicians did not consider suprachoroidal injection presenting new challenges, compared to other injections.**³
- No device / administration related complaints** were received across 4 recent clinical trials with three different drug products where current training program was implemented.⁷

Suprachoroidal injections are different from other ocular injection procedures, but specialized training improves physician confidence.

CONCLUSIONS & DISCUSSIONS

- Based on the criteria in Table 2 and the effectiveness assessment of the training (low complaint rates, increased physician comfort with procedure, and positive efficacy outcomes), this training program is rated at Level 4 for reducing possible rates of procedural complications, translating to **direct clinical benefits.**
- This **validated program is effective** in informing physicians of the critical tasks associated with suprachoroidal injections with the SCS Microinjector.

Clearside's proprietary training program, inclusive of the custom eye model and the training material content, has benefited from its wide use in investigational and commercial settings.

Clearside has refined the training process for suprachoroidal delivery using the SCS Microinjector® over the course of numerous clinical stage programs around the world.

References:
1. Henry CR, Walter SD, Chang PY, Warrow DJ, Naeini PE, Blinder KJ, Brevetti T, Yassine M, Dacey MS, Chu DS, Raiji VR, Rifkin LM, Shah M, Singer MA. Early adoption of triamcinolone acetonide suprachoroidal injection for uveitic macular edema: a physician survey. BMC Res Notes. 2024 Oct 23;17(1):317.
2. Nathan Fisher, Shelley Hancock, Thomas A Ciulla; Ophthalmic Procedure Training During COVID-19: Virtual and In-Person Training of the Suprachoroidal Injection Procedure. Invest. Ophthalmol. Vis. Sci. 2021;62(8):1190.
3. Wan CR, Kapik B, Wykoff CC, Henry CR, Barakat MR, Shah M, Andino RV, Ciulla TA. Clinical Characterization of Suprachoroidal Injection Procedure Utilizing a Microinjector across Three Retinal Disorders. Transl Vis Sci Technol. 2020 Oct 22;9(11):27.
4. FDA Guidance Document: Purpose and Content of Use-Related Risk Analyses for Drugs, Biological Products, and Combination Products, July 2024
5. Lee R, Raison N, Lau WY, Aydin A, Dasgupta P, Ahmed K, Haldar S. A systematic review of simulation-based training tools for technical and non-technical skills in ophthalmology. Eye (Lond). 2020 Oct;34(10):1737-1759.
6. McGaghie WC, Issenberg SB, Barsuk JH, Wayne DB. A critical review of simulation-based mastery learning with translational outcomes. Med Educ. 2014 Apr;48(4):375-85.
7. Clearside data on file, TRNG-0211

Related ARVO 2025 Posters:
4127-A0125, 592-A0229,
1780-B0515, 1782-B0517

Corresponding Author Contact:
cherry.wan@clearsidebio.com